

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922) 49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
Иваново (4932)77-34-06
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольяти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://kartell.nt-rt.ru> || kea@nt-rt.ru

НАКОНЕЧНИКИ К МИКРОПИПЕТКАМ



НАКОНЕЧНИКИ ОДНОРАЗОВЫЕ ДЛЯ АВТОМАТИЧЕСКИХ МИКРОПИПЕТОК

DESCRIPTION

Manufactured from virgin cadmium-free PP. High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum. The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered. These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market. Disposable Kartell tips come in three different package types to satisfy a variety of needs as economically as possible:

- bulk bag
- packed in racks
- refills for reloading racks

The indicated tip/micropipette compatibility may be subject to change, in case of any structural change carried out on micropipette models by manufacturers.




ПАКЕТ С ОДНОРАЗОВЫМИ ПИПЕТКАМИ



КОРОБКИ ДЛЯ ПИПЕТОК

DESCRIPTION

Some disposable tips are also available racked. Boxes containing the racks are in PP and have a snap-on lid. Boxes containing "macro" tips (5,000 μ l and 10,000 μ l) have a sliding horizontal cover. Containers and tips are autoclavable to +121°C for 20 minutes.



Система, состоящая из опорной плиты с множеством отверстий

DESCRIPTION

System comprising a support plate with an array of holes containing pipette tips. Designed for the reuse of empty racked tip boxes once the racked tips are used. It

can be used with different tip models, with the exception of the "macro" tips (5.000 µl and 10.000 µl). As an option for all refill tips, with the exception of "macro" tips (5.000 µl and 10.000 µl), spare empty racked tip boxes in PP are available.



НАКОНЕЧНИКИ ОБЪЕМ 0,1-10 МИКРОЛИТРОВ

DESCRIPTION

Compatibility: Gilson®, Eppendorf®, Socorex®, HTL®, Nichiryo®, Finnpiquette®, Thermo®

Manufactured from virgin cadmium-free PP.

High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Disposable Kartell tips come in different package types to satisfy a variety of needs as economically as possible:

- bulk bag
- packed in racks
- refills for reloading racks

PACKED IN RACKS

Some disposable tips are also available racked. Boxes containing the racks are in PP and have a snap-on lid.

Boxes containing "macro" tips (5,000 µl and 10,000 µl) have a sliding horizontal cover. Containers and tips are autoclavable to +121°C for 20 minutes.

REFILLS FOR RELOADING RACKS

System comprising a support plate with an array of holes containing pipette tips.

Designed for the reuse of empty racked tip boxes once the racked tips are used.

Several tip styles are available, with the exception of the "macro" tips (5.000 µl and 10.000 µl).

As an option for all refill tips, with the exception of "macro" tips (5.000 µl and 10.000 µl), spare empty racked tip boxes are available in yellow PP.

As an alternative, clear polycarbonate boxes with a horizontal sliding cover can be used; there are two different formats: one for refill tips with a volume range up to 250 µl and the other for refill tips with a volume range up to 1.000 µl.

Art.	Colour	Type	Packing
956	Neutral	Gilson® micro	Bulk

Art.	Colour	Type	Packing
9646	Neutral	Gilson® micro	Rack
9747	Neutral	Gilson® micro	Refill



НАКОНЕЧНИКИ ОБЪЕМ 0,5-20 МИКРОЛИТРОВ

DESCRIPTION

Compatibility: Kartell®, Eppendorf®, Gilson®, Biohit®, Brand®, Socorex®, Finnpiquette®, Thermo®, HTL®

Manufactured from virgin cadmium- free PP.

High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Disposable Kartell tips come in different package types to satisfy a variety of needs as economically as possible:

- bulk bag
- packed in racks
- refills for reloading racks

PACKED IN RACKS

Some disposable tips are also available racked. Boxes containing the racks are in PP and have a snap-on lid.

Boxes containing "macro" tips (5,000 µl and 10,000 µl) have a sliding horizontal cover. Containers and tips are autoclavable to +121°C for 20 minutes.

REFILLS FOR RELOADING RACKS

System comprising a support plate with an array of holes containing pipe tips.

Designed for the reuse of empty racked tip boxes once the racked tips are used.

Several tip styles are available, with the exception of the "macro" tips (5.000 µl and 10.000 µl).

As an option for all refill tips, with the exception of "macro" tips (5.000 µl and 10.000 µl), spare empty racked tip boxes are available in yellow PP.

As an alternative, clear polycarbonate boxes with a horizontal sliding cover can be used; there are two different formats: one for refill tips with a volume range up to 250 µl and the other for refill tips with a volume range up to 1.000 µl.

Art.	Colour	Type	Packing
959	Neutral	Eppendorf® Cristall	Bulk
9650	Neutral	Eppendorf® Cristall	Rack
9746	Neutral	Eppendorf® Cristall	Refill



НАКОНЕЧНИКИ ОБЪЕМ 2-20 МИКРОЛИТРОВ

DESCRIPTION

Compatibility: Eppendorf®, Gilson®, HTL®, Nichiryo®, Finnpiquette®, Thermo®

Manufactured from virgin cadmium- free PP.

High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Disposable Kartell tips come in different package types to satisfy a variety of needs as economically as possible:

- bulk bag
- packed in racks
- refills for reloading racks

PACKED IN RACKS

Some disposable tips are also available racked. Boxes containing the racks are in PP and have a snap-on lid.

Boxes containing "macro" tips (5,000 µl and 10,000 µl) have a sliding horizontal cover. Containers and tips are autoclavable to +121°C for 20 minutes.

REFILLS FOR RELOADING RACKS

System comprising a support plate with an array of holes containing pipe tips.


Designed for the reuse of empty racked tip boxes once the racked tips are used.

Several tip styles are available, with the exception of the "macro" tips (5.000 µl and 10.000 µl).

As an option for all refill tips, with the exception of "macro" tips (5.000 µl and 10.000 µl), spare empty racked tip boxes are available in yellow PP.

As an alternative, clear polycarbonate boxes with a horizontal sliding cover can be used; there are two different formats: one for refill tips with a volume range up to 250 µl and the other for refill tips with a volume range up to 1.000 µl.

Art.	Colour	Type	Packing
979	Neutral	Eppendorf®	Bulk
9657	Neutral	Eppendorf®	Rack
9745	Neutral	Eppendorf®	Refill



УНИВЕРСАЛЬНЫЕ ЖЕЛТЫЕ НАКОНЕЧНИКИ

DESCRIPTION

Capacity 2-200 µl

Compatibility: Gilson®, Eppendorf®, Socorex®, Nichiryo®, Biohit®, Brand®, Hamilton®
Kartell®

Capacity 5-200 µl

Compatibility: Kartell®, Eppendorf®, Gilson®, Biohit®, Brand®, Socorex®, HTL®, Nichiryo®,
Finnpipette®, Thermo®

Manufactured from virgin cadmium- free PP.

High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Disposable Kartell tips come in different package types to satisfy a variety of needs as economically as possible:

- bulk bag
- packed in racks
- refills for reloading racks

PACKED IN RACKS

Some disposable tips are also available racked. Boxes containing the racks are in PP and have a snap-on lid.

Boxes containing "macro" tips (5,000 µl and 10,000 µl) have a sliding horizontal cover. Containers and tips are autoclavable to +121°C for 20 minutes.

REFILLS FOR RELOADING RACKS

System comprising a support plate with an array of holes containing pipe tips.

Designed for the reuse of empty racked tip boxes once the racked tips are used.

Several tip styles are available, with the exception of the "macro" tips (5.000 µl and 10.000 µl).

As an option for all refill tips, with the exception of "macro" tips (5.000 µl and 10.000 µl), spare empty racked tip boxes are available in yellow PP.

As an alternative, clear polycarbonate boxes with a horizontal sliding cover can be used; there are two different formats: one for refill tips with a volume range up to 250 µl and the other for refill tips with a volume range up to 1.000 µl.

Art.	Colour	Type	Packing	Capacity
960	Yellow	Universal	Bulk	2-200 µl
9641	Yellow	Universal	Rack	5-200 µl
9749	Yellow	Universal	Refill	5-200 µl
980	Yellow grad.	Universal	Bulk	2-200 µl



НАКОНЕЧНИКИ ОБЪЕМ 5-200 МИКРОЛИТРОВ

DESCRIPTION

Compatibility: Gilson®, Kartell®

Manufactured from virgin cadmium-free PP.

High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Art.	Colour	Type	Packing
958	Yellow	Gilson®	Bulk



НАКОНЕЧНИКИ ОБЪЕМ 2-200 МИКРОЛИТРОВ

DESCRIPTION

Compatibility: Kartell®, Gilson®, Eppendorf®, Biohit®, Socorex®, HTL®, Nichiryo®, Finnpiquette®, Thermo®

Manufactured from virgin cadmium- free PP.

High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Disposable Kartell tips come in different package types to satisfy a variety of needs as economically as possible:

- bulk bag
- packed in racks
- refills for reloading racks

PACKED IN RACKS

Some disposable tips are also available racked. Boxes containing the racks are in PP and have a snap-on lid.

Boxes containing "macro" tips (5,000 µl and 10,000 µl) have a sliding horizontal cover. Containers and tips are autoclavable to +121°C for 20 minutes.

REFILLS FOR RELOADING RACKS

System comprising a support plate with an array of holes containing pipe tips.

Designed for the reuse of empty racked tip boxes once the racked tips are used.

Several tip styles are available, with the exception of the "macro" tips (5.000 µl and 10.000 µl).

As an option for all refill tips, with the exception of "macro" tips (5.000 µl and 10.000 µl), spare empty racked tip boxes are available in yellow PP.

As an alternative, clear polycarbonate boxes with a horizontal sliding cover can be used; there are two different formats: one for refill tips with a volume range up to 250 µl and the other for refill tips with a volume range up to 1.000 µl.

Art.	Colour	Type	Packing
940	Neutral grad.	Gilson®	Bulk
9643	Neutral grad.	Gilson®	Rack

Art.	Colour	Type	Packing
------	--------	------	---------

9751	Neutral grad.	Gilson®	Refill
------	---------------	---------	--------



НАКОНЕЧНИКИ ОБЪЕМ 5-200 МИКРОЛИТРОВ

DESCRIPTION

Compatibility: MLA®

Manufactured from virgin cadmium- free PP.


High clarity and excellent water-repellent properties designed to keep sample retention within the tip to a minimum.

The precise centering of the orifice at the end of the tip guarantees reliable directionality of the sample when delivered.

These tips are recommended for use with the Kartell automatic micropipettes but deliver excellent results even if used with the most other micropipette brands available on the market.

Art.	Colour	Type	Packing
------	--------	------	---------

968	Neutral	MLA®	Bulk
-----	---------	------	------



FINE TIP 0,1-200 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 2-300 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 20-300 МИКРОЛИТРОВ



УНИВЕРСАЛЬНЫЕ СИНИЕ НАКОНЕЧНИКИ



НАКОНЕЧНИКИ ОБЪЕМ 50-1000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 200-1000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 101-1000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 101-1000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 5-1000 МИКРОЛИТРОВ




НАКОНЕЧНИКИ ОБЪЕМ 5-1000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 1000-5000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 1000-5000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 500-5000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 2000-10000 МИКРОЛИТРОВ



ПУСТЫЕ КОРОБКИ ДЛЯ ПОПОЛНЕНИЯ

DESCRIPTION

Empty yellow PP boxes with snap-on covers for use with refill tips with capacity up to 1000 μ l.

MATERIAL

PP

Art.	Colour	Material	Type
9756-06	Yellow	PP	Standard lid



НАКОНЕЧНИКИ ОБЪЕМ 0,5-10 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 2-200 МИКРОЛИТРОВ



НАКОНЕЧНИКИ ОБЪЕМ 100-1000 МИКРОЛИТРОВ



СТЕРИЛЬНЫЕ ОДНОРАЗОВЫЕ НАКОНЕЧНИКИ ДЛЯ ПИПЕТОК С ФИЛЬТРОМ




НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 0,5-20 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 0,1-10 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 2-20 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 2-30 МИКРОЛИТРОВ




НАКОНЕЧНИКИ С ФИЛТРОМ ОБЪЕМ 2-100 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛТРОМ ОБЪЕМ 5-100 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛТРОМ ОБЪЕМ 2-200 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 20-300 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 50-1000 МИКРОЛИТРОВ



НАКОНЕЧНИКИ С ФИЛЬТРОМ ОБЪЕМ 50-1000 МИКРОЛИТРОВ



ПУСТЫЕ КОРОБКИ ДЛЯ НАКОНЕЧНИКОВ С ФИЛЬТРОМ ДЛЯ ПИПЕТОК

DESCRIPTION

Empty PP boxes with snap-on lid for pipette filter tips.

MATERIAL

PP

CERTIFICATION



Art.	Colour	Material
------	--------	----------

9756	Light Blue	PP
------	------------	----

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922) 49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
Иваново (4932)77-34-06
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольяти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://kartell.nt-rt.ru> || kea@nt-rt.ru